Optima Prone

User Manual

Be well, Live well
Wellell

CONTENTS

ENGLISH

USER MANUAL	1
-------------------	---

ESPAÑOL

MANUAL DE INSTRUCCIONES	26
-------------------------------	----

FRANÇAIS

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Model No.: 9P-086040

PLEASE READ ALL INSTRUCTIONS BEFORE USE

IMPORTANT SAFEGUARDS

READ ALL INSTRUCTIONS BEFORE USING

DANGER- To reduce the risk of electrocution:

- (1) Always unplug this product immediately after using.
- (2) Do not use while bathing.
- (3) Do not place or store this product where it can fall or be pulled into a tub or sink.
- (4) Do not place in or drop into water or other liquid.
- (5) Do not reach for a product that has fallen into water. Unplug immediately.

WARNING- To reduce the risk of burns, electrocution, fire, or injury to persons:

- (1) Evaluate patients for entrapment risk according to protocol and monitor patients appropriately.
- (2) The product may be used with caution for patients with spinal injury, but suggested to consult with physician before use. However, it should not be used for patients with unstable spinal fractures.
- (3) Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may result from a child swallowing a small part detached from the device.
- (4) Use this product only for its intended use as described in this manual. Do not use other mattress not recommended by the manufacturer.
- (5) Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to your supplier or Welld Inc. for examination and repair.
- (6) Keep the cord away from heated surfaces.
- (7) Never block any air openings of this product or place it on soft surfaces, such as a bed or couch, where openings may be blocked. Keep the air opening free of lint, hair, and other similar particles.
- (8) Never drop or insert any object into any opening or hose.
- (9) Do not modify this equipment without authorization of the manufacturer.
- (10) Mattress covers have passed skin sensitization and skin irritation test. However, If you suspect that you may have had or are having an allergic reaction, please consult a physician immediately.
- (11) Do not leave long lengths of tubing around the top of your bed. It could lead to strangulation.
- (12) Keep the pump away from flammable liquids or gases.
- (13) Do not use the extension cord for a long time.
- (14) Third conductor in the POWER SUPPLY CORD is only a functional earth.
- (15) The Operator cannot operate the instrument and touch the patient at the same time.

**CAUTION:**

If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (3.3m) between devices or turn off the mobile phone.

NOTE, CAUTION AND WARNING STATEMENTS:

NOTE: Indicate some tips.

CAUTION: Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.

WARNING: Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

SYMBOLS

	Authorized representative in the European community.		Consult operating instructions for use
	Catalog, reorder or reference number		Earth (ground)
	Manufacturer		Class II
	Complies with standards protecting against electric shock for type BF equipment.		Temperature limitation/temperature range
	Dry clean, Any Solvent Except Trichloroethylene		Do Not Bleach
	Do Not Iron		Do Not Dry Clean
	Tumble Dry, Normal, Low Heat		Machine wash, regular / normal, 95 degrees C (203 degrees F)
	Do Not Tumble Dry		Machine wash, regular / normal, 60 degrees C (140 degrees F)
	CE Marketing		Attention, should read the instructions.
	Protected against solid foreign objects of 12.5 mm and greater; Protection against vertically falling water drops		Medical Device
	Attention – Observe proper Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an appropriate collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.		

1. Introduction

This manual should be used for initial set up of the system and for reference purposes.

1.1 General Information

The system is a high quality mattress system suitable for treatment and prevention of pressure injury.

The system has been tested and successfully approved to the following standards:

E348970
53DG

MEDICAL EQUIPMENT- PUMP
WITH RESPECT TO ELECTRICAL SHOCK,
FIRE AND MECHANICAL HAZARDS ONLY
IN ACCORDANCE WITH ANSI/AAMI ES60601-1(2005). +
AMD(2012) and CAN/CSA C22.2 No. 60601-1(2014).

EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices to the IEC/EN 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.

1.2 Intended Use

This product is intended to help and reduce the incidence of pressure ulcers while optimizing patient's comfort. It can support patient up to 250 kg (550 lb). And it can also provide the following purposes:

- For long term home care of patients suffering from pressure ulcers.
- For pain management as prescribed by a physician.

The product can only be operated by personnel who are qualified to perform general nursing procedures and have received adequate training in knowledge of prevention and treatment of pressure ulcers.

NOTE: Equipment is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.

2. Patient Positioning Guide

Optima Prone is able to provide effective pressure relief to the patient being placed either in supine or in prone position.

WARNING

- (1) Please pay particular attention to the patient's head, neck and shoulders as each anatomical position should be placed on each corresponding section of mattress correctly. Upmost care must be provided to support the head and neck. Avoid any displacement or winding of tubes or lines.
- (2) Check regularly whether tubes and lines are in the correct position.
- (3) Obtain authorization from the clinician who is responsible for the patient's care before placing the patient in Prone Position.
- (4) When patient is being placed in the Prone position, areas that are most susceptible to pressure injuries, including forehead, chin, cheeks, ears, anterior shoulder, elbow, chest, genitalia, anterior pelvic bones (iliac crests & ischium), knees, dorsal feet and toes, should be treated with care¹.

1. European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed.). EPUAP/NPIAP/PPPIA: 2019. Section 8: Repositioning and Mobilization. Page 139 & 140.

2.1 General

Patient's shoulders should be aligned with the 4th cell, the shoulder support cell, either in the supine position or the prone position.

2.2 Supine Position

- (1) In Alternating Mode, all 21 cells will alternate according to the selected cycle time 10, 15, 20 or 25 minutes).
- (2) The 3 knobs in the Head Section are closed to support the head of the patient.

2.3 Activate the deflation of the head section when patient in supine position

- (1) The air cells in the head section can be deflated individually to provide nursing procedures such as hygiene, cannulation, etc. or provide pressure relief to patient's occiput.
- (2) This action should be conducted under the supervision of a competent clinician.
- (3) Make sure the head and neck are supported when deflating the air cells in the head section.
- (4) Under no circumstance should you leave the patient unattended while the air cells are deflated.
- (5) If the air cells in the head section are to stay deflated, the head and shoulders should be provided with appropriate support.

NOTE: The 4th cell of the mattress or the Shoulder Support cell will alternate when in alternating mode.

2.4 Prone Position– Not available for Homecare Environment

 **NOTE: Prone ventilation requires sufficient human resources and expertise to be performed safely.**

 **NOTE: Obtain authorization from the clinician who is accountable for the patient's care before placing the patient in Prone Position.**

- (1) Select Prone Position Mode to provide a firm and stable surface (Max Firm) for the following procedures.

 **NOTE: A stable surface is for the correct placement of the patient on the mattress.**

 **NOTE: Once the Prone Position Mode is selected, it will start the Prone Position timer automatically to show caregivers how long the patient has been positioned in prone position.**

- (2) Patient's shoulders should be aligned with the 4th cell, the Shoulder Support cell.
- (3) Support the head when turning the patient into the Prone Position.
- (4) Hyperextension of the neck should be avoided to provide a comfortable posture.

 **NOTE: A special design of the foam in the head section to relieve the pressure to the ear of patient.**

- (5) Avoid any displacement or winding of tubes or lines.
- (6) Make sure the patient's shoulders are still aligned with the 4th cell, the Shoulder Support cell.
- (7) Select a personalized repositioning schedule based on the professional assessment of the patient's condition by the competent clinician.

2.5 Reposition Patient's Head when in Prone Position – Not available for Homecare Environment

- (1) Select Head Reposition Mode to provide a firm and stable surface for the following procedures.
- (2) The air cell underneath the base cover will start inflating to lift the patient's shoulders and chest upward.
- (3) A beep sound indicates the air cell underneath the base cover is fully inflated.
- (4) Make sure the head and neck are supported when deflating the air cells in the head section.
- (5) Reposition the patient's head.
- (6) Turn the 3 knobs in the head section counter-clockwise to re-inflate the air cells.
- (7) Carefully place the patient's head on the fully inflated air cells.
- (8) Hyperextension of the neck should be avoided to provide a comfortable posture.

 **NOTE: A special design of the foam in the head section to relieve the pressure to the ear of patient.**

- (9) Avoid any displacement or winding of tubes or lines.
- (10) To cancel this function, select the Head Reposition mode again to cancel it and the air cell underneath the base cover will start to deflate slowly.

NOTE: Head Reposition Mode will automatically turn off after 10 minutes.


- (11) Make sure the patient's shoulders are still aligned with the 4th cell, the Shoulder Support cell.
- (12) The mattress will return to its previous therapy mode automatically to avoid long term pressurization of tissue.

3. Product Description

Unpack the box to check for any damage which may have occurred during shipment. If there are damages, please contact your dealer immediately.

3.1 Pump and Mattress System


- ❶ CPR
- ❷ Pump unit
- ❸ Mattress system
- ❹ Quick connector
- ❺ Position sensor

3.2 Pump Unit


Front

- ❶ Quick connector
- ❷ Position sensor
- ❸ Front panel


Rear

- ❶ Power switch
- ❷ Air filter
- ❸ Fuses
- ❹ Mounting brackets
- ❺ Battery cover

3.3 Front Panel

This model includes a supine position and a prone position. The front panel of the prone position is shown as follows:

(1) Power Key

Power Key is to turn on and turn off pump.

(2) Panel Lock key

User interface will be locked automatically 5 minutes after last operation. To unlock interface, press and hold Panel Lock key for 1 second.

⚠ NOTE: Always unlock user interface before turning off the pump.

(3) Max Firm key

This key is to make sure mattress can be reached maximal operating pressure. This mode will last for 20 minutes and then return to previous setting automatically. This function can be used during patient transfers or nursing procedures for better support. To disable this function, please press Max firm key again.

(4) LCD & Knob

Knob is to choose and confirm setting displayed on LCD, rotate the knob and press once to select the setting mode desired in the highlight.

(5) Alarms and alarm mute key

When there is an alert situation happened, the LCD screen will display a corresponding error code with an audible buzzer in an interval of a few seconds. The buzzer can be temporarily suspended by pressing the Alarm Mute key. Alarm buzzer will signal repeatedly every 3 minutes until the alert situation is resolved.

Below presents 2 types of alarm and please refer to the attached error code table for details.

3.3.1 For issues that can be checked and solved by user, such as low pressure alarm in head, body and foot areas, there'll be corresponding image. (see the table below)

(1) Low Pressure Alarm

Alarm will be indicated when the pump detects low pressure condition within the mattress.

- ⚠ **NOTE:** Under the alternating mode, for patient support, the system will first inflate and hold the cells with the right pressure level, keep monitoring the cells with low pressure condition. The system will remain in this state until the low pressure condition in the air cells is solved.
- ⚠ **NOTE:** The pump can seal the air inside the mattress, which temporarily keeps the patients from bottoming out during low pressure condition. Patient should be moved if the power failure issue cannot be solved soon.

(2) Low battery volume

- ⚠ **NOTE:** When low battery volume image is shown on the display, plug the power cord into the wall outlet immediately.
- ⚠ **CAUTION:** Pump system will stop running but display the alarm image and buzzer when low battery volume <10% is detected. To restore the system to the operating mode, plug the power cord into the wall outlet, and press the power key to restart the pump.

Code No.	Diagram	Description
01, 02, 05		Low pressure alarm: body section
03, 04, 06		Low pressure alarm: foot section
07		Low pressure alarm: head section
19		Position sensor is not connected

Code No.	Diagram		Description
20			Quick connector is not connected
21			Low battery volume (< 25 %)
21_1			Low battery volume (< 10 %)
21_2			Initial inflation: Low battery volume (< 10 %) After initial inflation: Low battery volume (≤ 2 %)
24			Connect a wrong mattress
30, 31, 34			Low pressure alarm: body section
32, 33, 35			Low pressure alarm: foot section
36			Low pressure alarm: Shoulder Lifting Mode

3.3.2 For issues that will affect the operation of pump and need to be solved by technician, it will show a wrench screwdriver with different error code. (see the table below)

(1) Power failure alarm


-  **NOTE:** There will be buzz sound and the alarm mute LED will light up constantly with dark screen when POWER FAILURE occurs.
-  **NOTE:** The pump can seal the air inside the mattress, which temporarily keeps the patients from bottoming out during power failure. Patients should be moved if the power failure issue cannot be solved soon.

Code No.	Diagram	Description
00		Initial inflation overtime
08-16		System malfunction
17		Pressure control system malfunction
18		System access control malfunction
37-41		System malfunction

4. Installation

4.1 Pump & Mattress Installation

- ⚠ **NOTE:** Unpack the box to inspect for any damage which may have occurred during shipment. If there are any damages, please contact your dealer immediately.
- ⚠ **NOTE:** USB function is provided for use by distributors or authorized technicians.

- (1) Place the mattress on top of bed frame and foot print toward foot end. Secure the mattress firmly by fixing the straps to bed frame to ensure mattress will be free from sliding. Securing straps are located on the base of the mattress.
 - ⚠ **NOTE:** Make sure that **ALL** the Relief knobs and emergency CPR knob are closed and locked in position.
 - ⚠ **CAUTION:** Placing another layers of sheet or bed linen between the patient and the mattress may reduce the benefits provided by the mattress and should be avoided at all time.
- (2) Hang the pump onto bed panel at the foot-end and adjust hangers to hang the pump to best upright position, or place the pump on a flat surface.
- (3) Connect the air hose connector of air mattress and the electrical cord of position sensor to the pump unit. Make sure the air hose connection is completed and secured when a "click" sound is felt or heard and position sensor assemble the alignment with the arrow pointing up.
 - ⚠ **NOTE:** Check and ensure the air hose is not kinked or tucked under mattress.
 - ⚠ **NOTE:** Auto recognition mattress of Optima Prone.

- (4)

⚠

NOTE: Make sure the pump unit is suitable for the local power voltage.

⚠

NOTE: The plug is also served to stop the device.

⚠

NOTE: Pump might still be in operation after the power cord disconnection if optional battery is built-in (9P-086520).
- (5)

⚠

CAUTION: The pump can only be applied to the mattress recommended by the manufacturer. Do not use it for any other purpose. (Applied part: air mattress) Do not position the equipment so that it is difficult to operate the disconnecting device.

⚠

CAUTION: Allow the mattress to inflate before positioning the patient on the mattress.
- (6)

The cable management tape is provided along the side of the mattress to manage power cord. Please make sure the wire has no signs of damage. After installation, extra length of the power cord, if any, should be neatly managed to prevent stumble accident and should be clear of moving bed mechanisms or other possible entrapment areas.

⚠

NOTE: The EQUIPMENT should be firmly placed at position where users can easily access to.

4.2 Position Sensor-installation & removal

Position sensor is key to auto seat inflation. It allows the mattress automatically to inflate when the head section has been raised to prevent bottoming out.
To initiate this function, you should connect the sensor cord to the unit.

To		Steps	Illustration
Attach	(1)	Align the arrows on the cord and the unit.	
	(2)	Insert the cord accordingly as illustrated.	
Detach	(1)	Put your fingers at the end of the cord as illustrated.	
	(2)	Pull backward to detach it from the unit. ⚠ Note: The grey area is the protection mechanism that prevents misuse and enables a proper detachment by sliding it backward.	

⚠

Warning: Do not immerse in water.

5. Operation


 **NOTE: Always read the operating instruction before use.**

5.1 General Operation

Switch the power switch on the side of the pump and press power key to turn it on.


5.1.1 Pre-set program

- (1) If the pre-set program is desired, you can choose TICK and press the knob, pump will automatically.


Perform pre-set
program:
Mode → Alternating
Cycle time → 10 mins.

- (2) If you have preferred setting, you can choose Mode and Cycle time to enter setting.

- (3) Once pump recognizes mattress size through quick connector, it will begin to inflate mattress which takes around 30 minutes to complete initial inflation.

- (4) After completing initial inflation, system will conduct automatic detection of patient weight and provide corresponding pressure to optimize mattress pressure relieving performance.

 **NOTE: Patient is NOT supposed to lie on the mattress before initial inflation is completed, but to be on the mattress when the auto adjustment process about to start.**

The 5 blue bars displayed on the screen represents the status of the automatic adjustment detection. When completes, the 5 bars will stop running and become


 **NOTE: Once the patient is on the mattress, the pump will automatically sense and adjust the pressure in the cells to support the patient.**


 **NOTE: Pressing the knob for a minimum of 3 seconds can quickly start the auto adjustment.**

- (5) Main LCD screen displays the current therapy. Setting can be changed by turning the


knob and the mode which is chosen will be highlighted, then press the knob once to enter or confirm the setting.


A. Mode

Alternating

Continuous Low Pressure

Back to previous page

NOTE: When changing between operating modes, the patient's monitoring and repositioning program should be reviewed.

B. Cycle time

10 mins

15 mins

20 mins

25 mins

Back to previous page


C. Tuning

Original setting

Decrease mattress pressure

Increase mattress pressure

Back to previous page


NOTE: The pressure in the air cells can be manually adjusted for patient comfort by entering the tuning mode in the lower right corner of the LCD screen.


NOTE: You can also check if the pressure is suitable for the patient by sliding one hand underneath the air cells at the level of the patient's buttocks. Always leave at least 2.54 cm (1 inch) space between patient and the Continuous Low Pressure cell to prevent bottoming-out.


D. Mode Choose

Supine Position Mode

Prone Position Mode

Back to previous page

E. Head Reposition mode


The bottom air cell will start to inflate to lift the patient's shoulder and chest upward for caregivers to reposition the patient's head easier.

F. Shoulder Lifting Mode

This mode can prevent hyperextension of the neck by lifting the patient's chest up to certain angle when in prone position.


Two angles are provided, medium and maximum, to fit each patient's condition.


No angle

Medium angle

Maximum angle


Back to previous page

G. Prone Position Timer

Once the Prone Position Mode is selected, the timer will start counting automatically.

5.2 Emergency CPR Operation

CPR knob is located at the patient's left-hand side of the mattress near the head area. When there is a need to perform CPR on the patient, quickly turn the CPR knob to release air from mattress. The quick connector on the pump unit can be disconnected for even faster deflation.

5.3 Relief Function


Turning the knob to deflate the cell enables the caregivers to float the patient's and take away the pressure from bony prominences.

- The 3 air cells in the head section are equipped with a knob.
- The 11 air cells in the torso section and the 6 air cells in the lower leg section cells are equipped with a knob.
- The single Shoulder Support (4th) cell has no knob and cannot be deflated.


5.3.1 For complex patient needs off-loading/pressure relief:

You can off-load more than one area on the patient for longer periods, but with the following restrictions:

- Select only two cells in the torso section to deflate.
- Select only one to two cell in the heel section to deflate.
- In the Prone position select three cells in the head section to deflate, or in the Supine position select three cells in the head section to deflate.

Closed

Open

Rotate the knob in a counterclockwise direction to the 12 o'clock position to deflate the cell.


NOTE: Deflating more than two adjacent cells may affect the support to the patient during the normal alternating cycle.

5.4 Evacuation Handles

NOTE: When an emergency occurs, the designated handles can be used for caregivers to quickly and safely evacuate the patient from the bed in the facility.



5.5 Transportation

NOTE: When there is a need of patient transportation, pump with the power cord can be disconnected from the mattress during transport. For optimum comfort, WELLELL recommends caregivers to change to Continuous Low Pressure mode and wait for 5 minutes before disconnecting. Cap on the quick connector as shown on the left to keep the air inside the mattress. The mattress remains inflated at least 24 hours depending on the setting. Optional battery model (9P-086520) can support at least 6 hours of treatment mode.

5.6 Notifications of parts replacement

When service message of battery/filter replacement (Code no. 22/23) shows up on the LCD screen, the device requires battery/filter replacement service. Please contact local representative or place of purchase for further support.

NOTE: To skip the message, please press the knob and confirm.

NOTE: Service message of battery/filter replacement which will only show up at the time the device is turned on.

5.7 Built-in batteries (Model number: 9P-086040)

For model with built-in batteries, it provides non-stop 6-hour therapy when in power outage. You can refer to top-right side of LCD display for battery status. Please read information below for battery operation.

- (1) Use AC power whenever possible. Connect the pump to AC power outlet and turn on the power switch to fully charge the battery for emergencies.
- (2) The warranty period is 6 months for the built-in batteries.
- (3) Connect the pump to AC power outlet immediately when the low-battery graphic shows up.
- (4) The battery may not be fully charged upon receipt. Connect the pump to AC power outlet for at least 10 hours.

 **NOTE: Pump cannot complete initial inflation when battery capacity is under 25%.**

6. Cleaning

It is important to follow the cleaning procedures to avoid cross contamination. Be sure to clean the surface in a dry and dust free environment. Wipe down the pump unit with a damp cloth pre-soaked with a mild detergent. Avoid contact with dust and proximity to dusty areas. Make sure that any cleaning agents you use will not harm or corrode the plastic casing on the pump unit. If your doctor or medical facilities have other special cleaning instruction, please follow the professional instruction.

 **CAUTION: Do not immerse or soak pump unit/position sensor.**

 **WARNING: Do not remove the casing of the pump to avoid electrical shock. All disassembly or repair should be done by authorized technicians.**

 **CAUTION: The pump does not need oil lubrication; please do not dissemble the system.**

Cover Material: • Stretch	    
Cover Material: • Nylon/PU	    

Wipe-down the mattress unit with a damp cloth pre-soaked with warm water containing a mild detergent, or chlorine bleach followed by an approved intermediate level disinfectant. Also the mattress top cover can be completely removed for laundry with water temperature indicated on the laundry label; however, it is recommended that the user still check with local policy to determine the time/temperature ratio required achieving thermal disinfection. The cover may also be cleaned using sodium hypochlorite diluted in water. After cleaning, please avoid dust and proximity to dusty areas and all parts should be air dried thoroughly before use.

 **CAUTION: Do not use phenolic based products for cleaning.**

 **CAUTION: Keep mattresses out of direct exposure to sunlight, and make sure to dry the mattress after cleaning.**

 **CAUTION: Make sure to execute the cleaning procedure every time before a new patient lying on the air mattress.**

7. Storage

- (1) To store the mattress, lay the mattress out flat and upsides down.
- (2) Roll from the head end towards the foot end with CPR valve open.
- (3) Foot-end strap can then be stretched around the rolled mattress to prevent unrolling.
- (4) After the product is stored for 3 months, it needs to be recharged again if optional battery is built-in (9P-086520).

NOTE: Do not kink, crease or stack the mattresses and do not store the system in direct sunlight, high temperature or moist areas.


8. Maintenance

8.1 General

- (1) Check main power cord and plug if there is abrasions or excessive wear.
- (2) Check mattress cover for signs of wear or damage. Ensure mattress cover and tubes are stubbed together correctly.
- (3) Check airflow from the quick connector. The airflow should alternate between two ventilation tubes in the middle and the other two at the bottom of quick connector every half-cycle time if it's in alternating mode.
- (4) Check the air hoses if there is kink or break. For replacement, please contact local representative or place of purchase.


8.2 Fuse Replacement

- (1) Disconnect the power plug when a blown fuse is suspected.
- (2) Remove the cover of the fuse holder by means of a flat screwdriver.
- (3) Insert a new fuse of the correct rating in, and replace the cover of the fuse holder back. The fuse should be rated as T1AL/250V type and VDE approved.

8.3 Air Filter Replacement

- (1) Open the air filter cap located at the back of pump.
- (2) Replace with a clean filter and put the cover lid back. The filter is reusable and can be washed gently with a mild detergent and water. Dry the filter before use.
- (3) Check and replace air filter regularly if environment is dirty. Suggested to do it every six month.

9. Expected Service Life

The products are intended to offer safe and reliable operation when use or installed according to the instructions provided by Wellell. Wellell recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function and indication on products. Otherwise, service and inspection of the devices generally should not be required.

10. Troubleshooting

Q1 Power is not ON

- Check if the plug is connected to mains.
- Check for a blown fuse.

Q2 Low Pressure Alarm is on

- Check if the Quick Connector is tightly secured.
- Check if all tubing connections along mattress are secured.
- Check if the air hoses are kinked or broken.
- Check if the CPR knob is tightly closed.

Q3 Power Failure Alarm is on

- Check if the power is suddenly shut down.
- Check if the power cord is connected properly.

Q4 Patient is bottoming out

- Pressure setting might be inadequate for the patient, increase comfort range 1 level higher and wait for a few minutes for best comfort.

Q5 Mattress form is loose

- Check if all the snap buttons or straps of mattress are all securely fastened.
- Check if the mattress is fixed to the bed frame by straps.

Q6 No air produced from some air outlets of the air tube connector

- This is normal since there is alternating mode. Air outlets take turns to produce air during their cycle time.

Q7 System malfunction (Code 11 or 12)

- Check if there is only one or two knobs at the foot-end of the mattress remained open.

If the above information does not solve your problems, please contact your local agent directly. They might require a technician to take care of the problem.

11. Technical Specifications

Item		Specification	
Power Supply (Note: See rating label on the product)		With battery	9P-086040: AC 110-120V 60 Hz, 1A
Fuse Rating		T1AL, 250V	
Cycle time		10, 15, 20, 25 minutes	
Dimension (L x W x H)		34.1 x 16.5 x 26.0 cm	
Weight		5 kg /11 lbs with battery	
Environment	Atmospheric Pressure	700 hPa to 1060 hPa	
	Temperature	Operation: 50°F to 104°F (10°C to 40°C) Storage: 5°F to 122°F (-15°C to 50°C) Shipping: 5°F to 158°F (-15°C to 70°C)	
	Humidity	Operation: 10% to 90% non-condensing Storage: 10% to 90% non-condensing Shipping: 10 % to 90% non-condensing	
Classification		Class II with functional earth, Type BF, IP21 Applied Part: Air Mattress ⚠ Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)	
Mattress		Specification	
Model		8" Mattress	
Dimension (L x W x H)		200 x 80/85/90 x 20 cm 78.7" x 31.5"/33.5"/35.4" x 8"	
Weight		14 kg / 30.8 lbs	
Max. Support Weight		250 kg / 550 lbs	
Relief Function		20 knobs(1 st to 3 rd , 5 th to 21 st air cells)	

NOTE:

- Consult the distributor or EU representative for further technical documents.
- The specifications are also suitable for other areas operating with same power supply.
- Mattress dimension and weight are measured without foam base.
- The manufacturer reserves the right to modify the specifications without notice.
- Cycle time may vary according to the local power supply.

Appendix A: EMC Information

Guidance and Manufacturer's Declaration- Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF emissions CISPR 11	Group1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	
 Warning: (1) The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used. (2) Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. (3) 3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the PUMP, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.		

Guidance and Manufacturer's Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Electrostatic Discharge (ESD) IEC61000-4-2	±8 kV contact ± 2, 4, 8, 15 kV air		±8 kV contact ± 2, 4, 8, 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC61000-4-4	±2 kV, 100 kHz for power supply line ±1 kV, 100 kHz for input/output line		±2 kV, 100 kHz for power supply line	Mains power quality should be that of atypical commercial or hospital

				environment
Surge IEC61000-4-5	± 0.5, 1 kV (line to line) ± 0.5, 1, 2 kV (line to earth)		± 0.5, 1 kV (line to line) ± 0.5, 1, 2 kV (line to earth)	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	Voltage Dips: i) 100% reduction for 0.5 period, ii) 100% reduction for 1 period, iii) 30% reduction for 25/30 period, Voltage Interruptions: 100% reduction for 250/300 period		120V	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	30 A/m	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF EM Fields IEC61000-4-3	3 V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz 380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM ±5 kHz deviation , (1kHz sine)	10 V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz 380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM ±5 kHz deviation, (1kHz sine)	3 V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz 380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM ±5 kHz deviation , (1kHz sine)	Recommended separation distance $d = \sqrt{P}$ 150kHz to 80MHz $d = 0.6\sqrt{P}$ 80MHz to 800MHz $d = 1.2\sqrt{P}$ 800 MHz to 2.7GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). ^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey , should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol:

				
NOTE 1: U _r is the a.c. mains voltage prior to the application of the test level				
NOTE 2: At 80 MHz and 800 MHz, the higher frequency range applies.				
NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people				
a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.				
b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.				

Recommended separation distances between portable and mobile RF communications equipment and this device:

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \sqrt{P}$	80 MHz to 800 MHz $d = 0.6\sqrt{P}$	800 MHz to 2,7 GHz $d = 1.2\sqrt{P}$
0.01	0.1	0.06	0.12
0.1	0.31	0.19	0.38
1	1	0.6	1.2
10	3.1	1.9	3.8
100	10	6	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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All rights reserved.
Wellell Inc.