

Sedens 500

Dynamic Seat Cushion

Lightweight seat cushion with long-lasting battery power

The Sedens 500 aligns with the 2019 clinical guideline, featuring advanced cushions and back supports designed to prevent pressure injury and manage pressure effectively. Lightweight and portable, it's perfect for chairs or wheelchairs, offering exceptional postural stability with a solid foam base and medial thigh support. Its long-lasting lithium-ion battery provides up to 12 hours of uninterrupted use, enhancing convenience. Visual and audible alarms prioritize user safety, making the Sedens 500 a comprehensive and innovative solution for mobility challenges while promoting comfort and well-being.



Optimal Pressure Redistribution

Effective alternating mode and static therapy mode to fulfill individual needs, and it also offers 5 comfort levels and 3 cycle times option for selection.



Extended Outdoor Use

The device enables extended outdoor use with its 12-hour lithium-ion battery on a full charge.

Sedens 500

Advanced mobility seat cushion: prolonged power supply and safety features



Specifications

Model	Sedens 500					
Pump	Dimension	9.8 x 4.7 x 2.2 in				
	Weight	2.2 lb				
	Supply voltage	100-240V / 50-60Hz				
	Therapy mode	Alternating / Static				
Cushion	Dimension	16" x 16" x 4"	17" x 17" x 4"	18" x 16" x 4"		
		18" x 18" x 4"	20" x 16" x 4"	20" x 18" x 4"		
	Type	Contour foam base with 4 air cells				
	Weight	2.2 - 2.6 lb				
	Top cover material	Poly / PU				
	Cell material	Nylon / PU				
	Base material	Nylon / PU				
	Maximum patient weight	264 lb				
HCPCS Code		E2624				
Ordering Information (For North America Only)						
	16" x 16"	17" x 17"	18" x 16"	18" x 18"	20" x 16"	20" x 18"
System Standard Cover	PD01905	PD01901	PD01902	PD01906	PD01903	PD01904
Pump	PD01001					

Pump: water resistant standards (IP2X); Mattress: flame retardant standards (EN597-1, EN597-2), RoHS, WEEE
Specifications are subject to change for improvement without notice

Wellell America Corp.
927 Mariner Street,
Brea, CA 92821, USA
T +1 714 671 3818
F +1 714 494 8173
info.usa@wellell.com
www.wellell.com
Print-2024/08 All rights reserved